

anti- A1BG antibody

Product Information

Catalog No.:	FNab10349
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Background	

Background

None

Immunogen information

Immunogen:	alpha-1-B glycoprotein
Synonyms:	A1B, A1BG, ABG, alpha 1 B glycoprotein, Alpha 1B glycoprotein, DKFZp686F0970, GAB, HYST2477
Observed MW:	40,70 kDa
Uniprot ID :	P04217


Application

Reactivity:	Human, Mouse			
Tested Application:	ELISA, WB, IHC			
Recommended dilution: WB: 1:500-1:2000; IHC: 1:50-1:200				
Image:				


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206) Tel :(0086)027-87384275





Immunohistochemistry of paraffin-embedded human kidney using FNab10349(A1BG antibody) at dilution of 1:50



Mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab10349(A1BG Antibody) at dilution of 1:1000

2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

www.fn-test.com