

anti- HK1 antibody

Product Information

Catalog No.:	FNab10353
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	hexokinase 1
Synonyms:	Brain form hexokinase, Hexokinase 1, Hexokinase type I, HK I, HK1, HK1 ta, HK1 tb, HK1 tc, HKI, HXK1
Observed MW:	100-120 kDa
Uniprot ID :	P19367

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:20 - 1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab10353(HK1 Antibody) at dilution of 1:100



Various lysates were subjected to SDS PAGE followed by western blot with FNab10353(HK1 Antibody) at dilution of 1:1000