

anti- ATP citrate lyase antibody

Product Information

Catalog No.:	FNab10354	
Size:	100µg	
Form:	liquid	
Purification:	Immunogen affinity purified	
Purity:	\geq 95% as determined by SDS-PAGE	
Host:	Rabbit	
Clonality:	polyclonal	
IsoType:	IgG	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

Immunogen information

Immunogen:	ATP citrate lyase
Synonyms:	ACL, ACLY, ATP citrate (pro S) lyase, ATP citrate lyase, ATP citrate synthase, ATPCL, Citrate cleavage enzyme, CLATP
Observed MW:	110-120kd
Uniprot ID :	P53396


Application

Reactivity:	Human, Mouse, Rat			
Tested Application:	ELISA, WB			
Recommended dilution: WB: 1:500 - 1:2000				
Image:				

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)





Various lysates were subjected to SDS PAGE followed by western blot with FNab10354(ATP citrate lyase Antibody) at dilution of 1:1000

2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275