

# anti- HK3 antibody

### **Product Information**

Catalog No.:	FNab10358	
Size:	100µg	
Form:	liquid	
Purification:	Immunogen affinity purified	
Purity:	$\geq$ 95% as determined by SDS-PAGE	
Host:	Rabbit	
Clonality:	polyclonal	
Clone ID:	None	
IsoType:	IgG	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

# Background

Hexokinases phosphorylate glucose to produce glucose-6-phosphate, the first step in most glucose metabolism pathways. This gene encodes hexokinase 3. Similar to hexokinases 1 and 2, this allosteric enzyme is inhibited by its product glucose-6-phosphate.

#### **Immunogen information**

Immunogen:	hexokinase 3 (white cell)
Synonyms:	Hexokinase 3, hexokinase 3 (white cell), Hexokinase type III, HK III, HK3, HKIII, HXK3
Observed MW:	None
Uniprot ID :	P52790

# Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution	:IHC: 1:50 - 1:200


#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com



# Image:



Immunohistochemistry of paraffin-embedded human colon cancer using FNab10358(HK3 antibody) at dilution of 1:100

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

www.fn-test.com