

anti- VAV3 antibody

Product Information

Catalog No.: FNab10378

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

This gene is a member of the VAV gene family. The VAV proteins are guanine nucleotide exchange factors (GEFs) for Rho family GTPases that activate pathways leading to actin cytoskeletal rearrangements and transcriptional alterations. This gene product acts as a GEF preferentially for RhoG, RhoA, and to a lesser extent, RAC1, and it associates maximally with the nucleotide-free states of these GTPases. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

Immunogen information

Immunogen: Guanine nucleotide exchange factor VAV3

Synonyms: VAV-3

Observed MW: 90-100 kDa Uniprot ID: Q9UKW4

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)


Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>

1



Recommended dilution: WB: 1:500-1:2000

Image:



mouse spleen were subjected to SDS PAGE followed by western blot with FNab10378(VAV3 antibody) at dilution of 1:800