

anti- CHM antibody

Product Information

Catalog No.: FNab10380

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen: Rab proteins geranylgeranyltransferase component A 1

Synonyms: DXS540, GGTA, HSD-32, REP-1, TCD, RAE1

Observed MW: 100 kDa, 60 kDa

Uniprot ID: P24386

Application

Reactivity: Human

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human colon cancer using FNab10380(CHM antibody) at dilution of 1:100



HEK293 cells were subjected to SDS PAGE followed by western blot with FNab10380(CHM antibody) at dilution of 1:500