

anti- CaMKII delta antibody

Product Information

Catalog No.: FNab10387

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen: calcium/calmodulin-dependent protein kinase II delta

Synonyms: CaM kinase II subunit delta, CaMK II subunit delta, CAMK2D,

CAMKD, CaMKII delta

Observed MW: 30-40 kDa Uniprot ID: Q13557

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:500; IF: 1:10-1:100

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human breast cancer using FNab10387(CAMK2D antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab10387(CAMK2D antibody) at dilution of 1:500