

anti- AIFM2 antibody

Product Information

Catalog No.:	FNab10395
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	apoptosis-inducing factor, mitochondrion-associated, 2
Synonyms:	AIFM2, AMID, Apoptosis inducing factor 2, ferroptosis suppressor protein 1, FSP1, p53 responsive gene 3 protein, PRG3, RP11 367H5.2
Observed MW:	41 kDa
Uniprot ID :	Q9BRQ8


Application

Reactivity:	Human		
Tested Application:	ELISA, WB		
Recommended dilution:WB: 1:200-1:2000			
Image:			

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)





2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

www.fn-test.com