

# anti- GREM1 antibody

### **Product Information**

Catalog No.: FNab10411

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## Immunogen information

Immunogen: Gremlin-1

Synonyms: C15DUPq, CKTSF1B1, CRAC1, CRCS4, DAND2, DRM, DUP15q,

GREMLIN, HMPS, HMPS1, IHG-2, MPSH, PIG2

Observed MW: 36 kDa
Uniprot ID: 060565

# **Application**

Reactivity: Human

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500 - 1:2000; IHC: 1:50 - 1:200

Image:


1

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human colon using FNab10411(GREM1 antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab10411(GREM1 antibody) at dilution of 1:1000