

# anti- BACE1 antibody

#### **Product Information**

Catalog No.:	FNab10412
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

## Background

Responsible for the proteolytic processing of the amyloid precursor protein(APP). Cleaves at the N-terminus of the A-beta peptide sequence, between residues 671 and 672 of APP, leads to the generation and extracellular release of beta-cleaved soluble APP, and a corresponding cellassociated C-terminal fragment which is later released by gamma-secretase.

## **Immunogen information**

Immunogen:	beta-site APP-cleaving enzyme 1
Synonyms:	BACE, KIAA1149
Observed MW:	50-55 kDa
Uniprot ID :	P56817

## Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution	:WB: 1:500 - 1:2000; IHC: 1:50 - 1:200


#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206) Tel :( 0086)027-87384275


www.fn-test.com



## Image:



Immunohistochemistry of paraffin-embedded rat brain using FNab10412(BACE1 antibody) at dilution of 1:100



Hela cells were subjected to SDS PAGE followed by western blot with FNab10412(BACE1 antibody) at dilution of 1:1000

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

www.fn-test.com