

anti- Cytokeratin 7 antibody

Product Information

Catalog No.: FNab10415

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen: keratin 7

Synonyms: CK 7, CK7, Cytokeratin 7, K2C7, K7, keratin 7, KRT7, Sarcolectin,

SCL, Type II keratin Kb7

Observed MW: 54 kDa Uniprot ID: P08729

Application

Reactivity: Human, Mouse

Tested Application: ELISA, IHC, WB

Recommended dilution: WB: 1:500 - 1:2000; IHC: 1:50 - 1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab10415(KRT7 Antibody) at dilution of 1:50



Hela cells were subjected to SDS PAGE followed by western blot with FNab02217(KRT7 antibody) at dilution of 1:1000