

# anti- CAB39L antibody

## **Product Information**

Catalog No.:	FNab10425
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

## Background

Component of a complex that binds and activates STK11/LKB1. In the complex, required to stabilize the interaction between CAB39/MO25(CAB39/MO25alpha or CAB39L/MO25beta) and STK11/LKB1(By similarity).

#### **Immunogen** information

Immunogen:	calcium binding protein 39-like
Synonyms:	Antigen MLAA 34, CAB39L, MO25 BETA, Mo25 like protein, MO2L, RP11 103J18.3
Observed MW:	40 kDa
Uniprot ID :	Q9H9S4

## Application


Reactivity:	Human, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution	:WB: 1:500 - 1:2000; IHC: 1:50 - 1:200

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206) Tel :( 0086)027-87384275



## Image:



Immunohistochemistry of paraffin-embedded human kidney using FNab10425(CAB39L antibody) at dilution of 1:100



rat kidney were subjected to SDS PAGE followed by western blot with FNab10425(CAB39L antibody) at dilution of 1:500

# 2

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

www.fn-test.com