

anti- FAK antibody

Product Information

Catalog No.:	FNab10444
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	1C4
IsoType:	IgG2c
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Focal adhesion kinase 1
Synonyms:	FADK, FADK 1, FAK, FAK1, Focal adhesion kinase 1, pp125FAK, Protein tyrosine kinase 2, PTK2, PTK2 protein tyrosine kinase 2
Observed MW:	100-120 kDa
Uniprot ID :	Q05397


Application

Reactivity:	Human, Mouse, Rat	
Tested Application:	ELISA, WB, IHC, IF	
Recommended dilution: WB: 1:500-1:3000; IHC: 1:20-1:200; IF: 1:20-1:200		
Image:		


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)





Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab10444(FAK Antibody) at dilution of 1:100



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab10444(FAK Antibody) at dilution of 1:1000.

2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

www.fn-test.com