

# anti- NFATC3 antibody

#### **Product Information**

Catalog No.: FNab10459

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

# **Background**

Acts as a regulator of transcriptional activation. Plays a role in the inducible expression of cytokine genes in T-cells, especially in the induction of the IL-2.

## **Immunogen information**

Immunogen: nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 3

Synonyms: NF AT4, NF ATc3, NFAT4, NFATC3, NFATX

Observed MW: 130-170 kDa

Uniprot ID: Q12968

# **Application**

Reactivity: Human

Tested Application: ELISA, IHC, WB

Recommended dilution: WB: 1:200-1:2000; IHC: 1:50-1:200

Image:


#### 1

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human kidney using FNab10459(NFATC3 antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab10459(NFATC3 Antibody) at dilution of 1:1000

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>