

---

## anti- 2019 nCOV S1 protein antibody

### Product Information

Catalog No.:	FNab10528
Size:	100µg
Form:	liquid
Purification:	protein A+G purified
Purity:	≥95% as determined by SDS-PAGE
Host:	mouse
Clonality:	monoclonal
Clone ID:	2A6
IsoType:	IgG
Storage:	PBS pH 7.4,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	S1-RBD
Uniprot ID :	P0DTC2

### Application

Reactivity:	2019 nCOV S protein
Tested Application:	ELISA, WB
Recommended dilution:	ELISA: Use at an assay dependent concentration(Suggestions: Use a concentration of 0.1 - 1 µg/ml); WB: 1:1000-1:10000

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)