

anti- SELENBP1 antibody

Product Information

Catalog No.:	FNab10530
Size:	100 μ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20 $^{\circ}$ C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Selenium-binding protein 1
Synonyms:	HEL-S-134P, LPSB, SBP56, SP56, hSBP
Observed MW:	55kDa
Uniprot ID :	Q13228

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human colon tissue slide using FNab10530(SELENBP1 Antibody) at dilution of 1:100



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab10530(SELENBP1 antibody) at dilution of 1:1000