
anti- CA19-9 antibody

Product Information

Catalog No.:	FNab10532
Size:	100µg
Form:	liquid
Purification:	protein A+G purified
Purity:	≥95% as determined by SDS-PAGE
Host:	mouse
Clonality:	monoclonal
Clone ID:	5G9
IsoType:	IgM
Storage:	PBS with 0.1% Proclin300,pH 7.4, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	purified Cancer Antigen 19-9
Synonyms:	CA19-9

Application

Reactivity:	Human
Application:	Recommended for sandwich immunoassays in ELISA,CLIA,GICA. Each laboratory should determine an optimum working titer for use in its particular application

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com