

anti- SLC7A11 antibody

Product Information

Catalog No.:	FNab10533	
Size:	100µg	
Form:	liquid	
Purification:	Immunogen affinity purified	
Purity:	\geq 95% as determined by SDS-PAGE	
Host:	Rabbit	
Clonality:	polyclonal	
IsoType:	IgG	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

Immunogen information

Immunogen:	solute carrier family 7, (cationic amino acid transporter, y+ system) member 11
Synonyms:	Amino acid transport system xc, CCBR1, Cystine/glutamate transporter, SLC7A11, Xct
Observed MW:	35kDa, 55kDa
Uniprot ID :	Q9UPY5

Application


Reactivity:	Human, Mouse, Rat			
Tested Application:	ELISA, WB			
Recommended dilution:WB: 1:500-1:2000				
Image:				

1

Wuhan Fine Biotech Co., Ltd. B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206) Tel :(0086)027-87384275

www.fn-test.com





2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275