
anti- 2019 nCOV N protein antibody

Product Information

Catalog No.:	FNab10601
Size:	1mg
Form:	liquid
Purification:	protein A purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	monoclonal
Clone ID:	5E8
IsoType:	IgG
Storage:	PBS buffer, pH 7.5,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	2019 nCOV N protein
Observed MW:	48 kDa

Application

Reactivity:	2019 nCOV
Tested Application:	ELISA, WB, Colloidal Gold-Based(detection antibody)
Recommended dilution:	ELISA: Use at an assay dependent concentration(Suggestions: Use a concentration of 0.1 - 1μg/ml)

Example for COVID-19 antigen Detection Kit (colloidal gold based), using finetest
FNab10600, FNab10601, P3003:

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com

