

anti- SIGLEC10 antibody

Product Information

Catalog No.:	FNab10618
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Sialic acid-binding Ig-like lectin 10
Synonyms:	Sialic acid-binding Ig-like lectin 10; SLG2; Siglec-10; Siglec-like protein 2; UNQ477/PRO940; sialic acid binding Ig like lectin 10
Observed MW:	76-100 kDa
Uniprot ID :	Q96LC7

Application

Reactivity:	Human, Rat
Tested Application:	WB
Recommended dilution:WB :	1:500-1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Jurkat cell lysates were subjected to SDS PAGE followed by western blot with FNab10618(SIGLEC10 antibody) at dilution of 1:1000