

## **anti- HBA1C antibody**

### **Product Information**

Catalog No.:	FNab10622
Size:	1mg
Form:	liquid
Purification:	Protein G purification
Purity:	≥90% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	8H2
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide, pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### **Immunogen information**

Immunogen:	hemoglobin A1c
Synonyms:	HbA1c, A1C, Hb1c

### **Application**

Reactivity:	Human
Tested Application:	particle-enhanced turbidimetric immunoassay

### **Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)