
anti- his-IF2DI tag antibody

Product Information

Catalog No.:	FNab10627
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.06% PC300 and 50% glycerol pH 7.3,-20℃ for 12 months(Avoid repeated freeze / thaw cycles.)

Background

This antibody can recognize recombinant protein with his-IF2DI tag.

Immunogen information

Immunogen: his-IF2DI protein

Application

Reactivity: his-IF2DI tag