

---

## anti- HBc antibody

### Product Information

|               |  |
|---------------|--|
| Catalog No.:  | FNab10631  |
| Size:         | 1mg  |
| Form:         | liquid   |
| Purification: | Protein A+G purification   |
| Purity:       | ≥95% as determined by SDS-PAGE   |
| Host:         | Mouse  |
| Clonality:    | monoclonal   |
| IsoType:      | IgG1   |
| Storage:      | PBS with 0.1% Proclin300 pH 7.3,store at 2-8° C up to 21 days, and at -20° C for long-term storage. (Avoid repeated freeze / thaw cycles.) |

### Immunogen information

Immunogen: hepatitis B virus core antigen

### Application

Reactivity: Recommended for ELISA, RIA and CLIA. Each laboratory should determine an optimum working titer for use in its particular platform application.