
anti- HBs antibody

Product Information

Catalog No.:	FNab10632
Size:	1mg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
IsoType:	IgG1
Storage:	PBS with 0.1% Proclin300 pH 7.3,store at 2-8° C up to 21 days, and at -20° C for long-term storage. (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen: hepatitis B virus surface antigen

Application

Reactivity: Recommended for ELISA, RIA and CLIA. Each laboratory should determine an optimum working titer for use in its particular platform application.

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com