

---

## anti- TLR9 antibody

### Product Information


Catalog No.:	FNab10635
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	Toll-like receptor 9
Synonyms:	Toll-like receptor 9, TLR9, CD289
Observed MW:	130 kDa
Uniprot ID :	Q9NR96

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded rat spleen tissue slide using FNab10635(TLR9 Antibody) at dilution of 1:100



HT-29 cells were subjected to SDS PAGE followed by western blot with FNab10635(TLR9 Antibody) at dilution of 1:1000