
anti- NEU3 antibody

Product Information


| | |
|---------------|---|
| Catalog No.: | FNab10636 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Immunogen information

| | |
|--------------|---|
| Immunogen: | N-acetyl-alpha-neuraminidase 3 |
| Synonyms: | NEU3, NEU3 protein, SIAL3, Ganglioside sialidase, Membrane sialidase, N-acetyl-alpha-neuraminidase 3, Sialidase-3 |
| Observed MW: | 48 kDa |
| Uniprot ID : | Q9UQ49 |

Application

| | |
|---------------------------|--------------|
| Reactivity: | Human |
| Tested Application: | ELISA, WB |
| Recommended dilution: WB: | 1:500-1:2000 |
| Image: | |



Various lysates were subjected to SDS PAGE followed by western blot with FNab10636(NEU3 Antibody) at dilution of 1:1000