
anti- IGF2 antibody

Product Information


Catalog No.:	FNab10641
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Insulin-like growth factor II
Synonyms:	Insulin-like growth factor II, Somatomedin-A, T3M-11-derived growth factor, IGF-II, IGF-2
Observed MW:	20 kDa
Uniprot ID :	P01344

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab10641(IGF2 Antibody) at dilution of 1:1000