

---

## anti- LIPG antibody

### Product Information


Catalog No.:	FNab10642
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	1G7
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Endothelial lipase
Synonyms:	LIPG, UNQ387/PRO719, Phospholipase A1, Endothelial cell-derived lipase, EDL, EL
Observed MW:	42-57 kDa
Uniprot ID :	Q9Y5X9

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10642(LIPG Antibody) at dilution of 1:1000