

---

## CDK5RAP2 antibody

### Product Information


Catalog No.:	FNab10643
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	CDK5 regulatory subunit-associated protein 2
Synonyms:	CDK5 regulatory subunit-associated protein 2, CDK5RAP2, CDK5 activator-binding protein C48, Centrosome-associated protein 215, CEP215
Observed MW:	250 kDa
Uniprot ID :	Q96SN8

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



Hela cells were subjected to SDS PAGE followed by western blot with FNab10643(CDK5RAP2 Antibody) at dilution of 1:1000