

FH antibody

Product Information

| | |
|---------------|---|
| Catalog No.: | FNab10646 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Protein A+G purification |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Mouse |
| Clonality: | monoclonal |
| Clone ID: | 1G6 |
| IsoType: | IgG2b |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Background


Also acts as a tumor suppressor.

Immunogen information

| | |
|--------------|---|
| Immunogen: | fumarate hydratase |
| Synonyms: | FH, Fumarase, fumarate hydratase, HLRCC, LRCC, MCL, MCUL1 |
| Observed MW: | 45-50 kDa |
| Uniprot ID : | P07954 |

Application

| | |
|--------------------------|-------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB |
| Recommended dilution:WB: | 1:1000-1:5000 |
| Image: | |



Various lysates were subjected to SDS PAGE followed by western blot with FNab10646(FH antibody) at dilution of 1:2000