
BTG3 antibody

Product Information


Catalog No.:	FNab10649
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Protein BTG3
Synonyms:	BTG3, ANA, TOB5, Abundant in neuroepithelium area protein, Protein Tob5, ANA, TOB5
Observed MW:	29 kDa
Uniprot ID :	Q14201

Application

Reactivity:	Human
Tested Application:	WB
Recommended dilution:WB:	1:500-1:1000
Image:	



Human liver were subjected to SDS PAGE followed by western blot with FNab10649(BTG3 Antibody) at dilution of 1:500