

HMOX1 antibody

Product Information


Catalog No.:	FNab10650
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	5G9
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	Heme oxygenase 1
Synonyms:	bK286B10, heme oxygenase (decycling) 1, Heme oxygenase 1, HMOX1, HO, HO 1, HO1, HO-1, HO-1/HMOX1, HSP32
Observed MW:	33 kDa
Uniprot ID :	P09601

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:100 - 1:200; IF: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab10650(HMOX1 antibody) at dilution of 1:100



Various lysates were subjected to SDS PAGE followed by western blot with FNab10650(HMOX1 Antibody) at dilution of 1:1000