
Catalase antibody

Product Information


Catalog No.:	FNab10653
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	9A6
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	catalase
Synonyms:	CAT, catalase
Observed MW:	60kd
Uniprot ID :	P04040

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB, IF
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab10653(Catalase antibody) at dilution of 1:100



Various lysates were subjected to SDS PAGE followed by western blot with FNab10653(Catalase antibody) at dilution of 1:2000