

---

## NLRC4 antibody

### Product Information


Catalog No.:	FNab10654
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	NLR family CARD domain-containing protein 4
Synonyms:	Caspase recruitment domain-containing protein 12, Ice protease-activating factor, NLRC4, CARD12, CLAN, CARD, LRR, and NACHT-containing protein (CED-4-like protein Clan, IPAF
Observed MW:	116 kDa
Uniprot ID :	Q9NPP4

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB, IF
Recommended dilution:	WB: 1:100-1:500; IF: 1:50-1:200
Image:	



Immunofluorescent analysis of THP-1 cells using FNab10654(NLRC4 antibody) at dilution of 1:50



N9 cells were subjected to SDS PAGE followed by western blot with FNab10654(NLRC4 antibody) at dilution of 1:400