
CST5 antibody

Product Information


Catalog No.:	FNab10655
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Cystatin-5
Synonyms:	Cystatin-D, CST5, CYTD
Observed MW:	16 kDa
Uniprot ID :	P28325

Application

Reactivity:	Human, Mouse
Tested Application:	WB
Recommended dilution: WB:	1:500-1:2000
Image:	



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab10655(CST5 antibody) at dilution of 1:500