

FABP6 antibody

Product Information


Catalog No.:	FNab10682
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	fatty acid binding protein 6
Synonyms:	FABP6, Fatty acid binding protein 6, Gastrotropin, GT, I 15P, I BABP, I BALB, I BAP, ILBP, ILBP3, Ileal lipid binding protein, ILLBP, Intestinal 15 kDa protein
Observed MW:	14 kDa
Uniprot ID :	P51161

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



A375 cells were subjected to SDS PAGE followed by western blot with FNab10682(FABP6 antibody) at dilution of 1:1000