

---

## Acetyl-Histone H3(K27) antibody

### Product Information


Catalog No.:	FNab10687
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	A synthetic peptide
Synonyms:	H3C1, H3C2, H3C3, H3C4, H3C6, H3C7, H3C8, H3C10, H3C11, H3C12
Observed MW:	15 kDa
Uniprot ID :	P68431

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:1000-1:5000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10687 at dilution of 1:2000