

Goat anti- monkey IgG (gamma)

Product Information

Catalog No.:	FNab10699
Size:	1 mg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Goat
Clonality:	polyclonal
IsoType:	IgG
Storage:	10mMPBS,pH 7.3, 4°C for 12 months.

Immunogen information

Immunogen:	monkey IgG gamma chain
------------	------------------------

Application

Reactivity:	Monkey
Tested Application:	ELISA, WB

Note: This product is for research use only, not for diagnostic or therapeutic use.