

anti- DAB2 antibody

Product Information

| | |
|---------------|---|
| Catalog No.: | FNab10718 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Immunogen information

| | |
|--------------|---|
| Immunogen: | disabled homolog 2, mitogen-responsive phosphoprotein(Drosophila) |
| Synonyms: | DAB2, Disabled homolog 2, DOC 2, DOC2, FLJ26626 |
| Observed MW: | 58 kDa, 96 kDa |
| Uniprot ID : | P98082 |

Application

| | |
|-----------------------|-----------------------------------|
| Reactivity: | Human |
| Tested Application: | ELISA, WB, IHC |
| Recommended dilution: | WB: 1:500-1:5000; IHC: 1:20-1:200 |
| Image: | |

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human prostate tissue slide using FNab10718(DAB2 antibody) at dilution of 1:100



HeLa cell lysates were subjected to SDS PAGE followed by western blot with FNab10718(DAB2 antibody) at dilution of 1:1000