

anti- ZNF746 antibody

Product Information

Catalog No.:	FNab10734
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Transcription repressor that specifically binds to the 5'-TATTTT[T/G]-3' consensus sequence on promoters and repress transcription of PGC-1-alpha(PPARGC1A), thereby playing a role in regulation of neuron death.

Immunogen information

Immunogen:	zinc finger protein 746
Synonyms:	FLJ31413, zinc finger protein 746, ZNF746
Observed MW:	50-55kDa,70-85 kDa
Uniprot ID :	Q6NUN9

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded rat testis tissue slide using FNab10734(ZNF746 antibody) at dilution of 1:100



HepG2 cell lysates were subjected to SDS PAGE followed by western blot with FNab10734(ZNF746 antibody) at dilution of 1:1000