

anti- EIF2C1 antibody

Product Information

Catalog No.:	FNab10769
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	eukaryotic translation initiation factor 2C, 1
Synonyms:	AGO1, Argonaute1, eIF 2C 1, EIF2C, eIF2C 1, EIF2C1, GERP95, hAgo1, Protein argonaute 1, Q99
Observed MW:	none
Uniprot ID :	Q9UL18

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, IF
Recommended dilution:	IHC: 1:20-1:200; IF: 1:10-1:100
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human colon using FNab0769(EIF2C1 antibody) at dilution of 1:50