

---

## anti- Perforin antibody

### Product Information

Catalog No.:	FNab10778
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	perforin 1(pore forming protein)
Synonyms:	Cytolysin, FLH2, HPLH2, P1, Perforin, Perforin 1, PFN1, PFP, PRF1
Observed MW:	
Uniprot ID :	P14222

### Application

Reactivity:	Human, Rat
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:20-1:200
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded rat lung tissue slide using FNab10778(perforin Antibody) at dilution of 1:100