
anti- LH antibody

Product Information

Catalog No.:	FNab10817
Size:	100µg
Form:	liquid
Purification:	protein A+G purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3H3
IsoType:	IgG1
Storage:	PBS with 0.1% proclin300, pH 7.4,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	purified luteinizing hormone
Synonyms:	Lutropin, Luteinizing hormone, LSH
Uniprot ID :	P01229

Application

Reactivity:	Human
Tested Application:	ELISA, CLIA, GICA (detection antibody)
Recommended dilution:	Each laboratory should determine an optimum working titer for use in its particular application.

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com