

anti- GATA3 antibody

Product Information

Catalog No.:	FNab10823
Size:	100µg
Form:	liquid
Purification:	protein A+G purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	2A5
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	GATA binding protein 3
Synonyms:	GATA binding factor 3, GATA binding protein 3, GATA3, HDR
Observed MW:	50 kDa
Uniprot ID :	P23771

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human breast cancer using FNab10823(GATA3 antibody) at dilution of 1:100



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab10823(GATA3 antibody) at dilution of 1:1000