

anti- FGF9 antibody

Product Information

Catalog No.:	FNab10839
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Fibroblast growth factor 9
Synonyms:	FGF9,FGF-9,GAF,HBFG-9,HBGF-9,SYNS3
Observed MW:	23 kDa
Uniprot ID :	P31371

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500 - 1:2000Image:


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab10839(FGF9 antibody) at dilution of 1:1000