

## anti- FLNC antibody

### Product Information

Catalog No.:	FNab10842
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	filamin C
Synonyms:	ABP 280, ABP 280 like protein, ABP L, ABP280A, ABPA, ABPL, Actin binding like protein, Filamin 2, Filamin C, FLJ10186, FLN C, FLN2, FLNC, Gamma filamin
Observed MW:	280-300 kDa
Uniprot ID :	Q14315

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Various lysates were subjected to SDS PAGE followed by western blot with FNab10842(FLNC antibody) at dilution of 1:1000