

anti- β HCG antibody

Product Information

Catalog No.:	FNab10847
Size:	100 μ g
Form:	liquid
Purification:	protein A+G purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	8G6
IsoType:	IgG1
Storage:	PBS with 0.1% proclin300, pH 7.4,-20 $^{\circ}$ C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	purified beta Chorionic gonadotropin
Synonyms:	Choriogonadotropin, CG, Chorionic gonadotropin
Observed MW:	None
Uniprot ID :	P0DN86

Application

Reactivity:	Human
Tested Application:	ELISA, CLIA, GICA (detection antibody)
Recommended dilution:	Each laboratory should determine an optimum working titer for use in its particular application.

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com