

anti- C1GALT1 antibody

Product Information

Catalog No.:	FNab10853
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	core 1 synthase, glycoprotein-N-acetylgalactosamine 3-beta-galactosyltransferase, 1
Synonyms:	B3Gal T8, Beta 1,3 galactosyltransferase, C1GALT, C1GALT1, Core 1 beta3 Gal T1, Core 1 O glycan T synthase, T synthase
Observed MW:	42 kDa
Uniprot ID :	Q9NS00

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:1000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



HepG2 cell lysates were subjected to SDS PAGE followed by western blot with FNab10853(C1GALT1 antibody) at dilution of 1:600