

## anti- PSM1 antibody

### Product Information

Catalog No.:	FNab10857
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Probably involved in the repair of mismatches in DNA.

### Immunogen information

Immunogen:	PMS1 postmeiotic segregation increased 1(S. cerevisiae)
Synonyms:	PMSL1
Observed MW:	60-100 kDa
Uniprot ID :	P54277

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



HL-60 cell lysates were subjected to SDS PAGE followed by western blot with FNab10857(PSM1 antibody) at dilution of 1:1000