

anti- CDKN1C antibody

Product Information

Catalog No.:	FNab10862
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	cyclin-dependent kinase inhibitor 1C
Synonyms:	BWCR, BWS, CDKN1C, KIP2, p57, p57Kip2, WBS
Observed MW:	57 kDa
Uniprot ID :	P49918

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded rat brain using FNab10862(CDKN1C antibody) at dilution of 1:50



mouse kidney lysates were subjected to SDS PAGE followed by western blot with FNab10862(CDKN1C antibody) at dilution of 1:1000